



Product Catalog

Premium Supplies for Hearing
Healthcare Professionals

Sanibel
Supply

www.sanibelsupply.com

Valid from September 2025



Welcome to Sanibel Supply's Product Catalog!

At Sanibel Supply, we blend the essence of sanitation and the clarity of audible care to support audiologists in their mission to provide pristine environments conducive to health and hearing. Our commitment is to ensure the highest standards of quality, enhancing patient well-being and auditory health.

In this catalog, you will find a comprehensive selection of our latest and most popular products, each carefully chosen to support the hearing patient journey. Whether you need hygiene-focused hearing aids, auditory support supplies, or health-centric audiologist tools, our range caters to every need with precision and care.

Explore our catalog and discover the exceptional quality of our offerings.

Enjoy your journey through our range!

The Sanibel Supply Team

Index

04	Probe Tips and Accessories
07	Disposable Ear Tips
12	Audiometry and Spareparts
19	Patient Preparation
21	Electrodes and Cables
26	Newborn Screening
28	Ear Wax Removal
32	Otoscopy and Otoscopy Supplies
34	Hearing Aid Fittings
36	Hearing Aid Maintenance and Care
43	Impression Materials and Supplies
47	Earplugs and Batteries
49	Balance Supplies
51	Infection Control
53	Thermal Printer Paper and Labels
55	List of abbreviations




Probe Tips and Accessories

Designed for professional audiologists, our high-quality probe tips are essential for accurate otoacoustic emission measurements. These tips are meticulously crafted to fit comfortably within the ear canal when used with ear tips, ensuring precise results by reducing external interferences. Our probe tips support personalized patient care and consistent diagnostic accuracy, reflecting our commitment to excellence.

Probe Tips and Accessories

Remember to order ear tips with you probe



(ABR) (OAE)

Sanibel Snap Probe Accessories

Designed for:

GSI: Novus.

Interacoustics: Sera.

MAICO: easyScreen.

Item	Description	Unit
8522759	Infant A	25 pcs
8522761	Infant C	25 pcs
8522758	Infant A	100 pcs
8522760	Infant C	100 pcs
8516436	Probe Waxguard	1 pc
8517889	Probe Waxguard	10 pcs
8522461	Assortment of accessories including Infant tip A and C and Waxguard	



(IMP)

Sanibel Probe Tip

Use with Sanibel ADI Ear Tips.

Designed for:

Amplivox: Otowaves.

GSI: Allegro.

Interacoustics: MT10.

Item	Color	Unit:
8002592	Transparent	1 pc



(IMP)

Sanibel Probe Tip

Use with Sanibel ADI Ear Tips.

Designed for:

Interacoustics: Titan.

MAICO: easyTymp.

Item	Color	Unit:
8001701	Transparent	1 pc



(ABR) (IMP) (OAE)

Sanibel Probe Tip Kit

Product includes: Probe Tip and Gasket.

Designed for:

GSI: Tymptstar Pro and Novus.

Interacoustics: AT235, AA222, Sera and Titan vers. 2.

MAICO: easyScreen and touchTymp.

Item	Unit:
8501570	3 pcs



(IMP)

Sanibel Probe Tip Adaptor Kit

Probe Tip, Washer, Damper.

No re-calibration required.

Use with Sanibel ADI Ear tips.

Designed for:

Interacoustics: AT235, AT235h, AA220, AA222 (blue)

MAICO: MI24, MI26, MI34 old series, MI44.

Item	Color	Unit:
8500430	Transparent/Black	3 pcs



(OAE)

Interacoustics Probe Tips ER10D


No re-calibration necessary.

Designed for:

Interacoustics: OtoRead.

MAICO: ERO-SCAN Pro.

Item	Color	Unit:
8012836	Black	4 pcs



(OAE)

Interacoustics Probe Tip Kit

Use with Sanibel ADI ear tips.

Designed for:

Interacoustics: Eclipse.

Item	Size	Color	Unit:
8012960	Ø3.3 mm	Transparent	8 pcs



(OAE)

MAICO Remote Probe Tip

Replacements for **MAICO** ERO-SCAN.

Item	Color	Unit:
8102340	White	4 pcs

The Importance of Cleaning Probe Tips

1. With clean probe tips, you prevent contamination from previous samples influencing your measurements.

2. Regular cleaning protects your probe tips from corrosion and blockages, enhancing their durability.

3. Removing biological and chemical contaminants safeguards both the user and the environment.

These points highlight the importance of regular cleaning to maintain high standards and safety in your work.



ABR IMP OAE

Sanibel Probe Cleaning Tool
Designed for:


Amplivox: Otowaves.

GSI: TympStar™ Pro, Novus, Allegro.

Interacoustics: Titan™/AT235™/AA222™/Sera™/MT10™.

MAICO: touchTymp MI22, MI23, MI24, MI26, MI34, MI36, MI44.

Item	Unit
8011351	1 pc




ABR IMP OAE

ProxySoft 3-in-1 Floss
Used to clean probe tips.

Size of Individual Floss Strand: 700 mm
Size of Plastic Threader: 76.2 mm
Size of FlexYarn Brush Section: 152.4 mm
Size of Soft Nylon Floss: 457,2 mm

Item	Unit
8003565	100 pcs



ABR IMP OAE

Thornton Bridge & Impant Floss
Used to clean probe tips.

Size of Individual Floss Strand: 254 mm
Size of Plastic Threader on Each End: 63.5 mm
Size of FlexYarn Brush Section: 127 mm.

Item	Unit
8003566	30 pcs

OAE


Sanibel Probe Tube Removal Tool Kit
Designed for:

GSI: Corti

Interacoustics: OtoRead.

MAICO: ERO•SCAN Plus.

Item	Description	Unit:
8104159	Probe Tubes	100 pcs
8513131	Tube Removal Tool	1 pc
8513132	Kit	100 pcs + 1 pc





Disposable Ear Tips

Our ear tips are indispensable tools in audiology, meticulously engineered to ensure precision and comfort for your patients. Crafted for a secure seal, they guarantee accurate test results by minimizing external interference. With a range of sizes and shapes tailored to diverse ear structures, we prioritize personalized care. Our commitment to versatility ensures consistent and reliable diagnostic outcomes.

6

Important Reasons to Use a New Ear Tip for Each Patient



By using new ear tips for each patient, you minimize the risk of cross-contamination. Despite thorough cleaning, it's challenging to completely decontaminate ear tips.



Using new ear tips minimizes the risk of poor sealing. Over time, ear tips degrade, making it harder to achieve a tight seal, ultimately leading to less accurate results.



Using new ear tips eliminates the risk of allergic reactions caused by residue from cleaning agents.



Simplify the process. Cleaning ear tips is time-consuming, so why not allocate your staff's time to more meaningful tasks?



Your patients will value it. They feel safer and more comfortable knowing that new ear tips are used.



Save time. Using a fresh ear tip ensures a quicker and more reliable seal.

Disposable Ear Tips



Sanibel ADI Ear Tips

For use with the following acoustic immittance instruments and otoacoustic emission systems with CLEAR probe tip.

Designed for:

Amplivox: All Otowave versions.

GSI: Novus, TympStar Pro, Corti and Allegro.

Interacoustics: Titan, Lyra, AT235, AA220, Eclipse, MT10, Sera, OtoRead II.

MAICO: easyTymp, touchTymp MI 22, MI 23, MI 24, MI 26, MI 34, MI 44, easyScreen, ERO•SCAN Plus.

Alternative to AKR (KR) / IA series: Old MAICO MI-series + RaceCar black probe versions together with adaptor kit (8500430).

Natus: Zodiac (not 901) and AccuScreen.

Item	Description	Unit
8511445	Various sizes for tympanometry and OAE	114 pcs
8501572	Various sizes for newborn screening	120 pcs



Sanibel Umbrella Ear Tips

For use with the following acoustic immittance instruments and otoacoustic emission systems with CLEAR probe tip.

Designed for:

Amplivox: all Otowave versions.


GSI: Novus, TympStar Pro, Corti and Allegro.

Interacoustics: Titan, Lyra, AT235, AA220, Eclipse, MT10, Sera, OtoRead II.

MAICO: easyTymp, touchTymp MI 22, MI 23, MI 24, MI 26, MI 34, MI 44, easyScreen, ERO•SCAN Plus. Alternative to AKR (KR). IA series: Old MAICO MI-series + RaceCar black probe versions together with adaptor kit (8500430).

Natus: Zodiac (not 901) and AccuScreen.

Item	Description	Color	Unit
8012985	7.5 mm	Green	25 pcs
8012987	9 mm	Red	25 pcs
8012989	10 mm	Blue	25 pcs
8012991	11 mm	Yellow	25 pcs
8012993	13 mm	Green	25 pcs
8012995	15 mm	Red	25 pcs
8012997	19 mm	Blue	25 pcs
8012999	22 mm	Yellow	25 pcs
8012986	7.5 mm	Green	100 pcs
8012988	9 mm	Red	100 pcs
8012990	10 mm	Blue	100 pcs
8012992	11 mm	Yellow	100 pcs
8012994	13 mm	Green	100 pcs
8012996	15 mm	Red	100 pcs
8012998	19 mm	Blue	100 pcs
8013000	22 mm	Yellow	100 pcs



Sanibel Mushroom Ear Tips

For use with the following acoustic immittance instruments and otoacoustic emission systems with CLEAR probe tip.

Designed for:

Amplivox: All Otowave versions.

GSI: Novus, TympStar Pro, Corti and Allegro.


Interacoustics: Titan, Lyra, AT235, AA220, Eclipse, MT10, Sera, OtoRead II.

MAICO: easyTymp, touchTymp MI 22, MI 23, MI 24, MI 26, MI 34, MI 44, easyScreen, ERO•SCAN Plus. Alternative to AKR (KR) IA series: Old MAICO MI-series + RaceCar black probe versions together with adaptor kit (8500430).

Natus: Zodiac (not 901) and AccuScreen.

Item	Description	Color	Unit
8013001	7 mm	Blue	25 pcs
8013003	8 mm	Yellow	25 pcs
8012969	9 mm	Green	25 pcs
8012971	10 mm	Red	25 pcs
8012973	11 mm	Blue	25 pcs
8012975	12 mm	Yellow	25 pcs
8012977	13 mm	Green	25 pcs
8012979	14 mm	Red	25 pcs
8012981	15 mm	Blue	25 pcs
8012983	19 mm	Yellow	25 pcs
8013002	7 mm	Blue	100 pcs
8013004	8 mm	Yellow	100 pcs
8012970	9 mm	Green	100 pcs
8012972	10 mm	Red	100 pcs
8012974	11 mm	Blue	100 pcs
8012976	12 mm	Yellow	100 pcs
8012978	13 mm	Green	100 pcs
8012980	14 mm	Red	100 pcs
8012982	15 mm	Blue	100 pcs
8012984	19 mm	Yellow	100 pcs

Disposable Ear Tips



Sanibel Flanged Infant Ear Tips

For use with the following acoustic immittance instruments and otoacoustic emission systems with CLEAR probe tip.

Designed for:

Amplivox: All Otowave versions.

GSI: Novus, TympStar Pro, Corti and Allegro.

Interacoustics: Titan, Lyra, AT235, AA220, Eclipse, MT10, Sera, OtoRead II.

MAICO: easyTymp, touchTymp MI 22, MI 23, MI 24, MI 26, MI 34, MI 44, easyScreen, ERO•SCAN Plus. Alternative to AKR (KR) / IA series: Old MAICO MI-series + RaceCar black probe versions together with adaptor kit (8500430).

Natus: Zodiac (not 901) and AccuScreen.

Item	Description	Color	Unit
8502837	Pediatric	Green	25 pcs
8012963	3-5 mm	Red	25 pcs
8012965	4-7 mm	Blue	25 pcs
8012967	5-8 mm	Yellow	25 pcs
8502838	Pediatric	Green	100 pcs
8012964	3-5 mm	Red	100 pcs
8012966	4-7 mm	Blue	100 pcs
8012968	5-8 mm	Yellow	100 pcs



Not available in EU

IA - Mushroom Ear Tips

For use with many acoustic immittance instruments and otoacoustic emission systems.

Designed for:

GSI: TympStar Pro, Allegro, Novus and Corti.

Interacoustics: (Clear probe): AT235, AA222, Titan, Otoread II and MT10.

Madsen/Otometrics: Zodiac (new), Accuscreen.

MAICO: touchTymp RaceCar, touchTymp MI24, MI26, MI34, MI36, easyTymp.

Path Medical: All Sentiero products.

Item	Description	Color	Unit
8520348	3-6 mm	Green	100 pcs
8520349	6 mm	Yellow	100 pcs
8520350	7 mm	Red	100 pcs
8520351	8 mm	Blue	100 pcs
8520352	9 mm	Green	100 pcs
8520353	10 mm	Yellow	100 pcs
8520354	11 mm	Red	100 pcs
8520355	12 mm	Blue	100 pcs
8520356	13 mm	Green	100 pcs
8520357	14 mm	Yellow	100 pcs
8520358	16 mm	Red	100 pcs
8520359	18.5 mm	Blue	100 pcs



Not available in EU

AKR - Mushroom Ear Tips

Use with the following acoustic immittance instruments and otoacoustic emission systems with black probe tip.

Designed for:

Amplivox: Otowave (all with white probe).

GSI: All tympanometers from GSI 20 - 39 + Tympstar and Allegro.

Interacoustics: AT235, A220, Eclipse, MT10 and Otoread (I).

Madsen/Otometrics: Capella 2, Zodiac 901 and Biologic AuDX Pro (+Flex).

MAICO: Old RaceCar, old MI22-36, MI 22, MI 23, MI 24, MI 26, MI 34 and MI 44.

Neurosoft: Neuro-Audio, Oscilla: TSM 300/400 + T800 series.

Item	Description	Color	Unit
8500570	3-6 mm	Yellow	100 pcs
8107510	4 mm	Red	100 pcs
8107511	5 mm	Blue	100 pcs
8500560	6 mm	Green	100 pcs
8500550	7 mm	Yellow	100 pcs
8500540	8 mm	Red	100 pcs
8500530	9 mm	Blue	100 pcs
8500520	10 mm	Green	100 pcs
8500510	11 mm	Yellow	100 pcs
8500500	12 mm	Red	100 pcs
8500490	13 mm	Blue	100 pcs
8107520	14 mm	Green	100 pcs
8500480	16 mm	Yellow	100 pcs



Sanibel Infant Ear Tips

For use with many acoustic immittance instruments and otoacoustic emission systems.

Designed for:

GSI: TympStar.

Interacoustics: OtoRead.

MAICO: ERO•SCAN.

Item	Description	Color	Unit
8013005	3 mm	Yellow	25 pcs
8002036	4 mm	Red	25 pcs
8013006	3 mm	Yellow	100 pcs
8002035	4 mm	Red	100 pcs



Not available in EU

IT - Select Series Ear Tips

Designed for:

Inventis: Flute, Clarinet, Viola and Timpani.

Resonance: R16m, R26m, R36m, R15c, R25c and R35c.

Item	Description	Color	Unit
8522402	8 mm	Green	100 pcs
8522403	10 mm	Red	100 pcs
8522404	12 mm	Blue	100 pcs
8522405	14 mm	Green	100 pcs
8522406	16 mm	Red	100 pcs

Disposable Ear Tips



Sanibel ANA – Ear Tips

Designed for:

Interacoustics: Eclipse.

Madsen/Otometrics: AccuScreen® Classic.

Natus: Echo-Screen®.

Item	Description	Color	Unit
8106363	4.0 mm	Blue	50 pcs
8106362	4.5 mm	Red	50 pcs
8106364	5.0 mm	Green	50 pcs
8106366	Infant tree	Yellow	50 pcs
8106365	Mushroom	Green	50 pcs



Sanibel AOT – Ear Tips

Designed for:

Madsen/Otometrics: AccuScreen® and Bio-logic® AuDX®.

Path Medical: All Sentiero products.

Item	Description	Color	Unit
8106367	3.7 mm	Yellow	50 pcs
8106368	4.0 mm	Light green	50 pcs
8106369	4.5 mm	Light pink	50 pcs
8106370	5.0 mm	Light blue	50 pcs
8106371	4-7 mm	Transparent	50 pcs



Not available in EU

Grason & Associates OT Select Ear Tips

OT Select Series Single Use Eartips were designed to fit two different variants of probe design.

Ear tips with two probe holes for use with **Otodynamics** ECP, UGS, SNS, SGS probes.

Item	Description	Color	Unit
8102632	3.5 mm	Blue	100 pcs
8102633	4.5 mm	Blue	100 pcs
8102634	5.5 mm	Blue	100 pcs
8102635	6.5 mm	Blue	100 pcs
8102636	4.8 mm	Blue	100 pcs
8102637	6 mm	Green	100 pcs
8102638	8 mm	Red	100 pcs
8102639	10 mm	Blue	100 pcs
8102640	12 mm	Green	100 pcs
8102641	14 mm	Red	100 pcs
8102642	16 mm	Blue	100 pcs



Not available in EU

Grason & Associates MO Select Ear Tips

Designed for:

MADSEN® OTOflex™.

Item	Description	Color	Unit
8102653	3-6 mm	Yellow	100 pcs
8102654	6 mm	Red	100 pcs
8102655	7 mm	Blue	100 pcs
8102656	8 mm	Green	100 pcs
8102657	9 mm	Yellow	100 pcs
8102658	10 mm	Red	100 pcs
8102659	11 mm	Blue	100 pcs
8102660	12 mm	Green	100 pcs
8102661	13 mm	Yellow	100 pcs
8102662	14 mm	Red	100 pcs
8102663	16 mm	Blue	100 pcs
8102664	18 mm	Green	100 pcs



Not available in EU

Grason & Associates OT Select Ear Tips

OT Select Series Single Use Eartips were designed to fit two different variants of probe design.

Ear tips with triangle holes for use with **Otodynamics** UGD, SND, SGD probes.

Item	Description	Color	Unit
8102644	4.8 mm	Blue	100 pcs
8102645	6 mm	Green	100 pcs
8102646	8 mm	Red	100 pcs
8102647	10 mm	Blue	100 pcs
8102648	12 mm	Green	100 pcs
8102649	14 mm	Red	100 pcs
8102650	16 mm	Blue	100 pcs



Sanibel OAE Screeners Ear Tips

Designed for:

Biologic: AuDx, Welch Allyn OAE screener.

Natus: Echo-Screen® III.

Item	Description	Unit
8515505	AuDX Premie 3 mm	20 pcs
8515506	AuDX Tree 3.5-5mm	20 pcs
8515507	AuDX Infant 6 mm	20 pcs
8515436	Adult OAE foam	20 pcs
8515435	Pediatric OAE foam	20 pcs
8515438	Jumbo OAE foam	20 pcs
8107885	AuDX Premie 3 mm	100 pcs
8107889	AuDX Tree 3.5-5mm	100 pcs
8107887	AuDX Infant 6 mm	100 pcs
8102681	Adult OAE foam	100 pcs
8102682	Pediatric OAE foam	100 pcs
8103123	Jumbo OAE foam	100 pcs



Audiometry and Spare parts

We offer various types of headsets that suit different testing methods and patient profiles. By combining our headsets with tailored consumables and spare parts, you can keep your headsets clean and comfortable for your patients.

Circumaural Headsets

(AUD)



RadioEar DD450 High Frequency Headset

The DD450 is a circumaural headset that can be used for audiometry testing in the standard audiometry range 125 Hz – 8 kHz as well as high-frequency testing up to 20 kHz. The DD450 complies with IEC-60645-1 and ANSI S3.6 up to 16 kHz. When using an IEC60318-1 coupler with a type 1 adaptor on an artificial ear, the standard allows for calibration using the RETSPLs.

RoHS compatible

Good ambient noise isolation

Gold-plated connectors

Capable of full range testing

Item	Description	Length
8520321	30°, 2×6.3 mm gold-plated mono jack	2.5 m
8520324	Straight, 2×6.3 mm gold-plated mono jack	2.5 m
8520325	90°, 1×3.5 mm gold-plated stereo mini jack	2.5 m
8520308	Straight, 1×3.5 mm gold-plated stereo mini jack	2.5 m




(AUD)



RadioEar DD65v2 Audiometric Headset

Item	Description	Length
8520327	Straight, 1×3.5 mm gold-plated stereo mini jack	2.5 m
8520328	30°, 2×6.3 mm gold-plated mono jack	2.5 m
8520330	Straight, 2×6.3 mm gold-plated mono jack	2.5 m
8520331	90°, 1×3.5 mm gold-plated stereo mini jack	2.5 m

(AUD)



RadioEar HDA200 Replacement Cable, Symmetric

For HDA200 and DD65 headphones.

Item	Description	Length
8004172	6.3 mm/ 30° jack plug	2 m

(AUD)



RadioEar Sanitation Kit

Replacement kit, containing set of new cushions and felt. Designed to maintain hygiene and prolong the life of DD65v2 and DD450 headsets.

Item	Unit
8524258	6 pcs

(AUD)



RadioEar Headband Cushion

For Sennheiser HDA200, RadioEar DD450, DD65v2 and Peltor audiometric headbands.

Item	Unit
8000770	1 pc

(AUD)



Ear Cushion Cover

Replaces disinfectants. Headphone covers for classic audiometric headphones.

Item	Description	Unit
8513404	Assorteret, 4.5"	100 pcs

Supra-aural Headsets & Supplies



(AUD)

RadioEar DD45 Adult Audiometric Headset

Supra-aural audiometric headset. DD45 audiometric headset 10 ohm, HBA headband.

Item	Description	Size	Length
8106340	30°, 2×6.3 mm gold-plated mono jack	Adult	2 m
8506471	30°, 2×6.3 mm gold-plated mono jack	Adult	5 m



(AUD)

RadioEar DD45 Contra Headset

Item	Description	Headband	Length
8010895	30°, 1×6.3 mm gold-plated mono jack	HBC	2 m
8105676	30°, 1×6.3 mm gold-plated mono jack	HB3045C	2 m



(AUD)

RadioEar DD45 Audiometric Headset

Item	Description	Size	Length
8512446	30°, 2×6.3 mm gold-plated mono jack	Pediatric	5 m
8105268	30°, 2×6.3 mm gold-plated mono jack	Adult	2 m
8511347	30°, 2×6.3 mm gold-plated mono jack	Adult	5 m



(AUD)

RadioEar DD45 Audiometric Headset

Headband: RE-7. Transducer: DD45.

Item	Description	Headband	Length
8010954	30°, 2×6.3 mm gold-plated mono jack	RE-7	2 m



(AUD)

RadioEar Foldable Audiometric Headband

Item	Description	Unit
8010882	Black	1 set



(AUD)

RadioEar DD45 Transducer

Audiometric earphones engineered to match the performance and function of the most popular earphones available for audiological testing. 10 ohm receiver. Easy installation. Reliable. RoHS compatible.

Item	Unit
8507268	1 pc

Supra-aural Headsets & Supplies



RadioEar RE-7 Headband

Item	Unit
8537014	1 set



RadioEar DD45/TDH Replacement Cable

Single, mono.

Item	Description	Length
8507205	6.3 mm/ 30° jack plug	2 m



RadioEar DD45™/TDH39™ Replacement Cable, Symmetric

Item	Description	Length
8507203	30° jack plug.	2 m
8507204	30° jack plug.	5 m



RadioEar Foam Pad for Contra Headset

Designed for e.g. TDH39c™ and DD45c™.


Item	Unit
8507292	10 pcs



RadioEar Headband Cushion

For RadioEar HBA foldable audiometric headband.

Item	Unit
8010886	1 pc



Ear Cushion Cover

Replaces disinfectants. Headphone covers for classic audiometric headphones.

Item	Description	Unit
8107419	Assorteret, 3"	500 pcs



RadioEar Cushion

For DD45™ and TDH39™ headsets.

Item	Description	Color	Unit
8507206	DD51	Black	1 pc
8603077	MX41	Black	1 pc



RadioEar TDH39/DD45 Earphone Cushion, model 51

For TDH39, 49, 50 & DD45 Earphones.

Item	Color	Unit
8507196	Red and Blue	1 set



Sanibel USB Cables A to B

Item	Length	Unit
8011241	2 m	1 pc
8011243	5 m	1 pc

Insert Earphone & Supplies



(AUD)

RadioEar IP30 Insert Earphone

The IP30 Insert earphone from RadioEar has been designed to meet the ISO 389-2 standard and matches the performance of the most popular insert earphones used in the industry today. These devices have proven to be quite accurate in a variety of testing conditions and are very effective in filtering out unwanted background noise. This makes them an ideal choice for non-ideal examination settings. Additionally the RadioEar IP30 Insert earphone does not require the use of bulky uncomfortable headbands.

The RadioEar IP30 Insert earphone is available in 10 ohm and 50 ohm impedances, as well as versions designed for use in AEP and contralateral applications.

Ideal for testing in non-ideal surroundings

Lightweight and compact

RoHS compatible

Item	Description	Ohm	Length
8101884	30°, 2x6.3 mm gold-plated mono jack	10 Ohm	2 m
8101891	30°, 2x6.3 mm gold-plated mono jack	50 Ohm	2 m
8101892	30°, 2x6.3 mm gold-plated mono jack	50 Ohm	2.9 m
8514834	30°, 2x6.3 mm gold-plated mono jack	10 Ohm	5 m
8514909	30°, 2x6.3 mm gold-plated mono jack	50 Ohm	5 m



(AUD)

RadioEar IP30 Contra Insert Earphone Cable

Plug type: 30°, 1x6.3 mm gold-plated mono jack.

Item	Length	Unit
8528812	2 m	1 pc



(AUD)

RadioEar IP30 Contra Insert Earphone

Item	Description	Length
8101894	Straight, 1x3.5 mm gold-plated mono mini jack	0.5 m
8108559	30°, 1x6.3 mm gold-plated mono jack	2 m



(ABR) (AUD)

Sanibel Foam Ear Tips

Designed for all insert phones.

Item	Description	Unit
8527912	Pediatric	50 pcs
8527913	Pediatric	100 pcs
8527915	Pediatric	750 pcs
8527918	Adult, medium	50 pcs
8527920	Adult, medium	100 pcs
8527921	Adult, medium	750 pcs
8535534	Adult, jumbo	24 pcs



(ABR) (AUD)

3M™ E-A-RLINK Foam Ear Tips

Designed for all insert phones.

Item	Description	Size	Unit
8500090	Adult/M	13.7 mm	50 pcs
8106226	Adult/M	13.7 mm	100 pcs
8500110	Pediatric/S	9.7 mm	50 pcs
8106227	Pediatric/S	9.7 mm	100 pcs
8500100	Jumbo/L	17.8 mm	24 pcs



(ABR)

3M™ E-A-RLINK Infant Ear Tips

Designed for all insert phones.

Item	Color	Size	Unit
8500300	Red/Black	3.5 mm	20 pcs
8500310	Red/Black	4.0 mm	20 pcs

Headsets and Supplies

(ABR) (AUD)

Sanibel Sound Tubes and Ear Tips Adaptor Kit

Fits all RadioEar IP30™ and E-A-RTone® versions. Includes Sanibel ear tip adaptors.

Item	Color	Unit
8100590	1 Blue and 1 Red	2 pcs

(ABR)

Sanibel Sound Tubes and Nipples Kit
Fits all RadioEar IP30™ and E-A-RTone® versions. Includes tube connectors.

Item	Color	Unit
8501072	1 Blue and 1 Red	2 pcs

(ABR) (AUD)

Tube Nipples

Fits all RadioEar IP30™ and E-A-RTone® versions. Use it with sound tubes 3A.

Item	Color	Unit
8500140	White	10 pcs

(ABR) (AUD)

Sanibel Ear Tip Adaptor

Fits all RadioEar IP30™ and E-A-RTone® versions.

Item	Color	Unit
8500580	Transparent	6 pcs

(ABR)

Sanibel IP30 Tube Connectors

Fits all RadioEar IP30™ and E-A-RTone® versions.

Item	Color	Unit
8107391	Black	10 pcs

(ABR)

Sanibel IP30 Clip

Replacement velcro clip for insert earphones. Designed for a variety of insert earphones, including: IP30, E-A-RTone 3A and 5A, Etymotic ER-3A and ER-3C.

Item	Description	Unit
8537295	Clip	1 pc
8507410	Kit Includes velcro neckstrap and clip Set	

(ABR) (AUD)

Sanibel Sound Tubes

Single/Contra. Fits all RadioEar IP30™ and E-A-RTone® versions. The product include adaptor nipples.

Item	Unit
8107396	1 pc

(ABR) (AUD)

Sanibel Ear Tip Adaptor Cleaning Tool

Small brush for cleaning the Sanibel Probe Tip Adaptors.

Item	Unit
8500830	10 pcs

Bone Transducer Headsets & Spare parts



RadioEar B71 Bone Transducer Headset
10 ohm.

Item	Description	Length
8506488	Straight, 1x6.3 mm gold-plated mono jack	2 m
8519447	3.5 mm plug	2 m
8510370	30°, 1x6.3mm gold-plated mono jack	2 m
8511359	30°, 1x6.3mm gold-plated mono jack Kit with 10 pcs. Headband	2 m



RadioEar B-81 Bone Transducer Headset
12.5 ohm.

Item	Description	Length
8509840	Straight, 1x6.3mm gold-plated mono jack	2 m
8510371	30°, 1x6.3mm gold-plated mono jack	2 m
8511361	30°, 1x6.3mm gold-plated mono jack Kit with 10 pcs. Headband	2 m



RadioEar Replacement Cable for Bone conductor B71/B81
With safety fixation point.

Item	Unit
8515713	1 pc



RadioEar Replacement Cable for old version of B-71
6.3 mm/ 30° jack plug.

Item	Length	Unit
8508473	2 m cable	1 pc




RadioEar Replacement Pads (3288-1)
For RadioEar Spring Steel Headband.

Item	Unit
8504777	10 pcs



RadioEar AMBAND Bone Conductor Headband (P1188)
For adult, children and infants.

Item	Unit
8507214	10 pcs




RadioEar Spring Steel Headband (P3333)

Item	Unit
8507213	1 pc

RadioEar Covered Headset for B71 & B-81 P3333C

Item	Unit
8506712	1 pc



Patient Preparation

As an audiologist, you play a vital role in ensuring accurate hearing measurements for your patients. Your expertise and care are crucial for a successful procedure. By selecting the appropriate products for patient preparation, including electrode gel, skin prep gel, and tape, you can ensure reliable results. Your dedication to patient comfort and well-being creates a positive experience and effective treatment.

Patient Preparation



(ABR) (OAE)

Softa® Swabs

Alcohol Pads, The Softa Swabs from B. Braun are ready-to-use swabs for local skin disinfection. 70% Isopropyl alcohol.

Item	Size	Unit
8500050	30×30 mm	100 pcs



(ABR) (OAE)

Gauze SoftaSwabs

100% pure cotton, highly absorbent and soft. Perfect for hygiene. Non-sterile, designed for non-invasive procedures..

Item	Size	Unit
8500040	75×75 mm	100 pcs



(ABR)

SignaGel® Electrode Gel

Bacteriostatic conductive saline electrode gel. Water soluble, non-staining, non-gritty.

Item	Unit
8102628	250 g



(ABR)

Spectra 360® Electrode Gel

Salt- & chloride-free conductive electrode gel suitable for long-term applications. Product of choice for conductive rubber/carbon electrodes.

Item	Unit
8500740	250 g



(ABR)

Everi Abrasive Paste

Reduces the skin impedance during EEG, EP, ECG and Holter Stress test examination.

Item	Unit
8012952	160 g



(ABR)

Sonaville® Electrode Gel

Formaldehyde-free and glycerol-free.

Item	Unit
8500380	250 ml



(OAE) (ABR)

Q-tips Cotton Tipped Applicators

Cotton tipped applicators for patient preparation. Non-sterile.

Item	Unit
8500810	1000 pcs



Paper Kidney Tray

Item	Unit
8520139	1 pc



(ABR) (OAE)

Ambu® Skin Prep Pads 2121M

An abrasive pad ideal for skin preparation, is designed for easy application on the finger due to its adhesive properties and shape.

Item	Size	Unit
8523414	19 × 28 mm	500 pcs



(ABR)

3M™ Red Dot™ Trace Prep

One-step skin preparation. Designed to remove non-conductive skin layer to achieve better signal quality.

Item	Size	Unit
8500030	18 mm × 5 m	1 pc



(ABR)

3M™ Transpore™ Surgical Tape

Hypoallergenic and not made with natural rubber latex. Water-resistant and porous for breathability. Easy to handle with gloves. Transparent.

Item	Size	Unit
8500750	25 mm × 9 m	12 pcs



(ABR)

3M™ Micropore™ Surgical Tape

Hypoallergenic and not made with natural rubber latex. Comfortable and allows the skin to breathe.

Item	Size	Unit
8500760	25 mm × 9 m	12 pcs



(ABR)

NuPrep™ Skin Prep Gel

Specialty gel with mild abrasive to clean and reduce the impedance of the skin before applying the electrodes.

Item	Description	Unit
8500320	114 g	3 pcs



(ABR)

Ten20™ Conductive EEG Paste

Adhesive gel for disc electrodes.

Item	Description	Unit
8500330	3×8 oz/228 g	3 jars



(ABR)

Ten20™ Conductive EEG Paste

Adhesive gel for disc electrodes.

Item	Description	Unit
8012924	3×4 oz/114 g	3 tubes



Electrodes and Cables

Our electrodes are indispensable in audiological tests like ABR and ECoG, capturing precise brain or muscle activity. Featuring hydrogel for optimal conductivity, our electrodes are available in a variety of sizes and brands to suit different clinical needs. They enable accurate data collection critical for diagnosing auditory and vestibular conditions. For dependable results, ensure proper placement and preparation.

Enhance your electrode measurements with our meticulously designed electrode cables, essential for transferring responses to your evoked potentials device. Explore our comprehensive range on the next pages to ensure reliable performance in your diagnostics.

Electrodes and Cables



(ABR)

Sanibel Snap Electrodes

Pre-gelled with a very gentle hydrogel, perfect for ABR testing on infants. Packed in 20 smaller pouches of 3 electrodes. Intended use areas: Newborn, diagnostic follow-up, ABR, C-Vemp.

Item	Size	Unit
8107137	2.2×3 cm	60 pcs




(ABR)

Sanibel Tab Electrodes

Pre-gelled cloth electrode with a very gentle hydrogel, perfect for ABR testing on infants. Packed in 20 smaller pouches of 3 electrodes. Intended use areas: Newborn, diagnostic follow-up, ABR, C-Vemp.

Item	Size	Unit
8107410	2.2×3 cm	60 pcs




(ABR)

Sanibel TM Electrode

For ECoG. Use with cable 8500650.

Item	Unit
8011489	10 pcs
8049550	1 pc



Choose the correct electrode for optimal performance

Sanibel Supply® offers various kinds of surface electrodes for ECG as well as EEG and ABR. Surface electrodes vary in size, adhesion, connectivity and frequency (test category).

The frequency plot above shows the electrical stimuli sent directly through the system compared to the surface electrode (here ABR). The two lines are similar and demonstrate that the surface electrode can be used for ABR testing.

At Sanibel Supply®, we always test and validate our surface electrodes to ensure optimal performance for your testing!




(ABR)

Sanibel Cable for TM ECoG Electrode

For Eclipse™ and Bio-logic® AEP system. 1.5 mm/.060" female connector. For use with Sanibel TM electrodes.

Item	Unit
8500650	1 pc



(ABR)

Sanibel ECoG Starter Kit

For Eclipse™ and Bio-logic® AEP system. Includes cable, gel and 2 pcs TM electrodes.

Item	Unit
8013209	4 pcs



(ABR)

Sanibel Gold Ear Tips Electrode

For ECoG. Use with cable 8011343.

Item	Description	Unit
8500440	Pediatric, 10 mm	20 pcs
8500370	Adult, 13 mm	20 pcs



(ABR)

Spes Medica Tip Trode Cables

For use with Sanibel Gold Ear Tips Electrodes. 1.5 mm with touch proof connectors and plated clamp clips.

Item	Length	Unit
8011343	50 cm	2 pcs

Electrodes and Cables



ABR

Sanibel Pinch Clip Cable

For use with snap electrodes.
1.5 mm/.060" female connector.

Item	Color	Length	Unit
8500390	RWB	50 cm	3 pcs
8500400	WYK	50 cm	3 pcs
8519048	WYK	100 cm	3 pcs
8500410	RWBYK	50 cm	5 pcs

ABR

Sanibel Alligator Clip Cable

For use with tab electrodes. 1.5 mm/.060" female connector.

Item	Color	Length	Unit
8535775	RWB	50 cm	3 pcs
8535776	WYK	50 cm	3 pcs
8535774	RWBYK	50 cm	5 pcs



ABR

Sanibel Universal Electrode Cable

Non-tangling.

Item	Color	Length	Unit
8011353	RWB	50 cm	3 pcs
8011344	RWBYK	50 cm	5 pcs

Not available in EU

ABR

Sanibel Button Electrode Cable

For use with snap electrodes.
Non-tangling.

Item	Color	Length	Unit
8532948	RWB	50 cm	3 pcs
8532949	RWBYK	50 cm	5 pcs

ABR

Spes Medica Disc Electrode

For use with Ten20™ electrode paste.

Item	Color	Length	Unit
8011498	RWB	50 cm	3 pcs
8011345	RWBYK	50 cm	5 pcs

ABR

Gold-Plated Cup Electrodes

Ø 10 mm gold-plated.

Item	Color	Length	Unit
8107145	RWBYK	50 cm	5 pcs



ABR


Sanibel EP Electrode Jumper Cable

Allows two of the preamplifier inputs to be connected, e.g. for linking mastoid electrodes for cortical testing.

Item	Unit
8107872	1 pc



Electrodes and Cables



Vermed Tab Electrode

Rectangular tab electrode. Conductive Adhesive Hydrogel.

Item	Size	Unit
8516225	12.7×22.86 mm	100 pcs



ABR

Vermed Wetgel Snap Electrode

Pre-gelled hydrogel electrode with Ag/AgCl sensor. Very suitable for diagnostic follow-up. Packaged 3 pcs. per pouch.

Item	Unit
8516221	Box/75 pcs



ABR

Ambu® Neuroline 720

The Ambu Neuroline 720 electrodes, which are pre-gelled and self-adhesive, are perfect for sleep studies. Their robust yet skin-friendly adhesion guarantees that the electrode remains secure throughout the entire examination, even in the presence of patient perspiration.

Item	Size	Unit
8500060	45 ×22 mm (non-wire)	25 pcs


Electrodes and Cables



Vermed Economy Line Snap Wetgel Electrode

Non-woven cloth electrode with conductive wet gel.

Item	Size	Unit
8516251	1.57x1.4"	50 pcs



Vermed Hydrogel Snap Electrode

Large soft electrode size making it very suitable for newborn screening and ABR measurements.

Item	Size	Unit
8516907	22.22x34.92 mm	60 pcs



Vermed Hydrogel Snap Electrode

Smaller than the rectangular version, but otherwise similar. Good, soft electrode size making it very suitable for newborn screening and ABR measurements.

Item	Size	Unit
8102022	0.875 in/22.22 mm	60 pcs




Vermed Economy Line Solid Gel Snap Electrode

Solid gel type snap electrode suitable for diagnostic follow-up. Latex free. Intended use areas: Newborn, diagnostic follow-up all ages.

Item	Unit
8516252	50 pcs





Newborn Screening

We offer specialized newborn hearing screening supplies so you can provide precise screening and optimal care for your youngest patients. Browse our range of screening consumables on the next page.

Newborn Screening

ABR

Sanibel Infant EarCups Kits

Designed for:

GSI: Novus.

Interacoustics: Titan, Sera.

Madsen/Otometrics: Accuscreen.

MAICO: easyScreen, MB 11 Classic.

Use only with equipment calibrated for over-ear style disposables. Disposable over-the-ear earphones are high quality and ideal for bilateral newborn hearing screenings. 1 set consists of 2 Infant EarCups™, 3 snap or tab electrodes and 1 alcohol pad. For use with a range of different providers of newborn hearing screening equipment.

Item	Description	Unit
8049600	Tab	20 sets
8049700	Snap	20 sets
8508848	w/o electrodes	20 sets

Sanibel Infant EarCup w/ Hydrogel

New Hydrogel version to ensure a more comfortable solution for the baby while testing. Choose any electrode type you wish to add to meet your preference. Order electrodes separately.

Product includes 1 set consists of 2 Infant EarCups™, 3 snap or tab electrodes and 1 alcohol pad.

Item	Unit
8516095	20 sets



ABR

Sanibel Infant EarCup Adaptor

Designed for:

GSI: Novus.

Interacoustics: Titan, Sera.

Madsen/Otometrics: Accuscreen.

MAICO: easyScreen, MB 11 Classic.

Connects insert phones with Infant EarCups™. NB: Do not contain filters.

Item	Description	Unit
8500450	EarCup Adaptor	2 pcs
8100577	Kit Includes adaptors and Sound tubes	2 pcs

ABR

Sanibel EARTurtle Adapter for RadioEar IP30

Item	Unit
8515872	1 pc





Ear Wax Removal

Ear wax and debris-free ear canals are necessary for reliable hearing testing. Browse our range of tools for ear wax removal and ear cleaning.

Ear Wax Removal

Bionix Safe Ear Curettes

Bionix® Safe Ear Curettes allow healthcare providers to safely and gently remove cerumen with styles to match different types of wax and personal preference. Bionix® Safe Ear Curettes are flexible, to reduce injury to the ear canal, making them safer than stainless steel curettes. Seven color coded tip styles to meet specific needs, giving you “the right tool for any job” for every cerumen situation.

Item	Description	Color	Size	Unit
8515585	MicroLoop®	Green	3 mm	50 pcs
8515586	ControLoop®	Orange	4 mm	50 pcs
8515587	AngleLoop®	Red	4 mm	50 pcs
8519359	VersaLoop®	Purple	3 mm	50 pcs
8519360	InfantScoop®	Blue	2 mm	50 pcs
8519372	CeraSpoon®	Yellow	4 mm	50 pcs
8519361	FlexLoop®	White	4 mm	50 pcs
8519362	Variety Kit, (10 each style & 15 white)			



Why does ear wax accumulate?

Disruption of the self-cleaning mechanism of the ear causes cerumen to accumulate, causing different problems within our ears. These can be related to e.g. general health of the ear as well as hearing.

Anatomical variations, infectious or dermatological diseases can cause narrowing or obstruction of the ear canal. This can then interfere with the normal migratory process of foreign agents from the ear.

Indications for treatment and treatment options

When symptomatic, cerumen accumulation requires treatment. Usually, ear wax is defined as impacted, if it disrupts visualization

of the tympanic membrane, as well as being connected with symptoms.

Often, removal of the cerumen is also done to allow a needed view of the tympanic membrane, e.g. for diagnostic purposes and audiometry.

The first-line option for treatment of symptomatic cerumen in patients is cerumenolytics (softening agents, e.g. eardrops). The use of cerumenolytic agents, such as eardrops, increase the probability of cerumen clearance.

Patients with heavier impaction or ear canal disease may require irrigation or manual removal under microscope by a trained doctor or an ENT specialist.



Bionix® Lighted Ear Curettes

Bionix Lighted Ear Curettes combine illumination and magnification to help healthcare providers take less time to see and safely remove cerumen from the ear canal. The Lighted Ear Curette with magnification brilliantly illuminates the ear canal, eliminating the practice of blindly curettaging the sensitive area. Complete visualization shortens procedure time and increases efficacy.

All versions include 1 Light Source and 1 magnifying lens.

Item	Description	Size	Unit
8519364	Lighted AngleLoop®	4 mm diameter	50 pcs
8519365	Lighted MicroLoop®	3 mm diameter	50 pcs
8519366	Lighted WaveCurette™	2.5 mm wide tip	50 pcs
8519367	Lighted FlexLoop®	4 mm wide tip	50 pcs
8519368	Lighted CeraSpoon®	4 mm wide tip	50 pcs
8519369	Lighted InfantScoop®	2 mm wide tip	50 pcs
8519373	Lighted VersaLoop®	3 mm wide tip	50 pcs
8519370	Lighted CeraPik™, Great for hard, impacted cerumen and foreign bodies	2.5 mm concave tip	50 pcs



All versions include 1 Light Source and 1 Magnifying Lens

Ear Wax Removal

Bionix OtoClear® Spray Wash Kit

A disposable ear irrigation system that is simple to use, effective for compacted cerumen removal and eliminates messy clean-up.


Item	Description	Unit
8519371	Kit of 20× irrigator tips, 1× 470ml spray wash bottle, 1× ear basin	
8519390	OtoClear® Ear Basin	3 pcs
8519375	Bionix Otoclear Irrigation Tips	40 pcs



Bionix Lighted Ear Curette Variety Pack

Bionix® Lighted Ear Curettes combine illumination and magnification to help healthcare providers take less time to see and safely remove cerumen from the ear canal. The Lighted Ear Curette with magnification brilliantly illuminates the ear canal, eliminating the practice of blindly curettaging the sensitive area. Complete visualization shortens procedure time and increases efficacy.

Item	Description
8519363	With 1 light source and magnifying lens. Includes 8 pcs. of the Lighted Flexloop and 7 pcs. of each of the other styles (not including Cerapik).



Bionix Lighted Forceps

Bionix Lighted Forceps combine illumination and magnification to help healthcare providers take less time to see and safely remove ear canal and nasal passage obstructions.

Item	Description
8519392	10 forceps, 1 light and 1 magnifying lens

Introducing EarWay® Pro

– your go-to solution for effective, safe, and hassle-free earwax removal.

This innovative tool ensures a quiet and noiseless procedure that takes less than a minute per ear. EarWay Pro allows for seamless ear wax management and continuation of care within the same appointment. Experience the convenience and efficiency that EarWay Pro brings to ear care.

Recommended by leading audiologists across the globe, the EarWay Pro gives autonomy back to clinicians.



EarWay Pro – Wax Removal Helix

Very easy and safe way to remove ear wax. The helix grabs around the ear wax lump and pulls it out in a comfortable way.

Item	Color	Size	Unit
8518024	Blue	4.2 mm	25 pcs
8518025	Purple	5.4 mm	25 pcs
8521854	Green	6.5 mm	25 pcs



Otoscopes and Otoscopy Supplies

Our range of Otoscopy supplies includes disposable ear specula and otoscope tips, ensuring that each examination is safe, hygienic, and maintains the functionality of our equipment.

Otoscopes and Otoscopy Supplies



(AUD) (HAS) (IMP)

Luxamed Otoscope Auris LED 2.5V

Innovative otoscope in pocket size. Automatic switchoff after 3 minutes. Life expectancy of LEDs: approx. 100.000 hours by constant current drive. incl. 20 disposable ear funnels, grey, each 10×2.5 mm and 10×4.0 mm. Designed for **Heine/KaWe** speculums.

Item	Unit
8519433	1 pc



Luxamed Otoscope Auris LED 3.7V - USB/PC Charged

Innovative otoscope in pocket size. Automatic switchoff after 3 minutes. Life expectancy of LEDs: approx. 100.000 hours by constant current drive incl. recharger with USB plug (adapter EU/UK/US). incl. 20 disposable ear funnels, grey, each 10×2.5 mm and 10×4.0 mm. Designed for **Heine/KaWe** speculums.

Item	Unit
8519434	1 pc



(AUD) (HAS) (IMP)

HEINE Mini 3000 LED Otoscope

Experience HEINE quality in a pocket-sized format. This compact otoscope offers exceptionally bright illumination and true-to-life color rendering with LEDHQ technology, eliminating the need for LED replacement. Its exclusive design and sturdy housing ensure durability and convenience, making it always close at hand. Available in black or blue, this otoscope combines intuitive operability with a detailed view and flexible power source options.

Item	Unit
8603028	1 pc



(AUD) (HAS) (IMP)

HEINE Mini-e Ear Light LED

This durable pocket-sized examination light provides a powerful, concentrated beam thanks to its unique LEDHQ technology. Ergonomically designed for comfort, it automatically switches off when placed in your pocket. With 20,000 guaranteed on/off cycles and replaceable AAA batteries, this light is both practical and reliable.

Item	Unit
8603029	1 pc



(HAS)

EARNimals

Cartoon faces that are placed on the otoscope. Make ear examinations more fun for kids.

Item	Description	Unit
8011501	Assorted	5 pcs



(HAS)

Specula Wall Dispenser

For ear funnels. Can contain 6 different sizes.

Item	Unit
8515678	1 pc



(AUD) (HAS) (IMP)

1 - Otoscope Specula

Disposable ear speculums.

Item	Description	Size	Unit
8507594	Economy for Heine/KaWe	2.5 mm	100 pcs
8507618	Economy for Heine/KaWe	4 mm	100 pcs
8514368	For Heine/KaWe	2.5 mm	100 pcs
8514367	For Heine/KaWe	4 mm	100 pcs
8514370	Genuine for Heine/KaWe	2.5 mm	50 pcs
8514369	Genuine for Heine/KaWe	4 mm	50 pcs
8107750	For Welch Allyn	2.5 mm	250 pcs
8108252	For Welch Allyn	4 mm	250 pcs

(AUD) (BAL) (HAS) (IMP)

2 - Sanibel Disposable Speculas for Video Otoscope

Designed for the **AURICAL** OTocam 300, **Interacoustics** VIOT and **XION** video otoscopes.

Item	Unit
8500470	50 pcs

(AUD) (BAL) (HAS) (IMP)

3 - Video Otoscopes Specula

Designed for the **Interacoustics** Viot and **OrlVision** otoscope.

Item	Description	Size	Unit
8515083	For KaWe	2.5 mm	100 pcs
8511728	For KaWe	4 mm	100 pcs



Heine® Otoscope Hard Case

Protect your otoscope from damage with a durable, hard shell case. For Heine Mini 2000 or 3000 Otoscopes.

Item	Unit
8536592	1 pc



Heine® Bulb

For Heine mini-c earlight


Item	Unit
8530892	1 pc



Replacement Bulb

For Heine mini3000 Otoscope.

Item	Unit
8515448	1 pc



Hearing Aid Fitting

Our REM tubes are integral to Real-Ear Measurement (REM), facilitating the precise fitting and verification of hearing aids. These high-quality tubes enable precise adjustments tailored to the unique hearing profiles and ear canal acoustics of each patient.

Hearing Aid Fitting



Sanibel REM Probe Tubes

Universal REM Probe Tubes.

Universal fit, length 76 mm, beveled 45° angle tip and black marker ring.

Designed for:

AudioScan

Auditdata

Interacoustics



Inventis

MedRx

Natus Aurical

Use: Rem Tubes.

Item	Description	Color	Unit
8533301	Type A	Black	50 pcs
8533303	Type C	Blue	50 pcs
8533304	Type D	Red	50 pcs



Sanibel Silicone Tubes for SPL60

Designed for:

Interacoustics: Affinity™ REM probe.

Use: Rem Tubes.

Item	Color	Unit:
8500700	Transparent	50 pcs



Sanibel Aidapter for Hearing Aid Fitting

Clear Silicone Acoustics couplers. A unique silicone piece that clicks onto the HA-ladapter. You can use the Aidapter with receiver-in-canal (RIC) and thin tube hearing aids. Compared to acoustic putty, the Aidapter™ provides a good acoustic seal and standardizes the coupling, ensuring repeatable coupler measurements. The Aidapter™ also eliminates the hassle and hygienic issues of using acoustic putty.


Item	Description	Unit
8525568	Interacoustics	3 pcs
8525578	MedRx	3 pcs



Sanibel™ REM Probe Tubes, type A

Universal fit, length 76 mm, beveled 45° angle tip and black marker ring.

Item	Unit
8104116	50 pcs



Ear Coupler Sealing Wax

Use to create a seal between coupler and hearing instrument.

Item	Unit
8500080	1 pc



Battery Adaptor

Designed for:


Interacoustics Affinity™, Affinity™ Compact, Callisto™, Equinox™.

MedRx: Avant Hit+.

Signia: UNITY 2.

Item	Description	Color	Unit
8011198	BAA675	Blue	1 pc
8011199	BAA13	Orange	1 pc
8011200	BAA312	Brown	1 pc
8011201	BAA10	Yellow	1 pc
8011202	BAA5	Red	1 pc





Hearing Aid Maintenance and Care

Proper maintenance and customization of customers' hearing aids are of paramount importance. We offer a wide range of products, including cleaning tools, disinfectants, spare parts, and accessories, all designed to facilitate this process. With access to these products, you can ensure that the hearing aids function optimally and are tailored to customers' individual needs and preferences.

Hearing Aid Maintenance and Care




(HAS)

Westone LED Earlite

The 18 Lumens light is concentrated in a small tip that fits most ears. Includes 2 AAA batteries and 2 earlite tips.

Item	Unit
8602609	1 pc



(HAS)

Westone Tube Puller

Allowing quick and simple tube replacement.

Item	Description	Unit
8527877	Monofilament	1 pc
8527878	Wire loop	1 pc



(HAS)

Westone Vent Cleaning Tool

Item	Unit
8527910	1 pc



(HAS)

Westone Cerumen Removal Tool

Cerumen Removal Tool features a loop on one end, a brush on the other and a magnet in the handle to help remove and installing hearing instrument batteries.

Item	Description	Unit
8527907	w/Brush and Magnet	1 pc
8527908	w/Brush	1 pc
8527909	Removal Tool	1 pc



(HAS)

Luxamed LED Ear Light With Acrylic Scaled Top

Bright LED light for placement of ear dams. Uses (2) AAA batteries (included). Black.

Item	Description	Unit
8519357	Luxamed LED Ear light with acrylic scaled top	1 pc
8519358	Acrylic Top With Scale	1 pc



(HAS)

Drewe Otoplastik Wire Cleaning Brush

Clean hearing instruments in an easy way.

Item	Unit
8519100	1 pc



Luxamed Stethoclip

Hearing Instrument listening stethoscope for reassuring the fitting process

Item	Description	Unit
8515600	Stethoclip, Black	1 pc
8520346	Ear funnels	1 pc
8520347	ITE adaptor	1 pc



Drewe Otoplastik Stethoclip With Coupling

Item	Description	Unit
8530900	ITE adaptor with tubing	1 pc

Hearing Aid Maintenance and Care



Westone Battery Extractor Tool


The hearing aid battery extractor tool unique design assists in removing even the smallest batteries from the tightest battery compartments. Quantity pricing discounts are available.

Item	Unit
8527911	1 pc



Westone Reamer for 13 Tube


Item	Unit
8536590	1 pc



Westone Probe Set

For probing, material removal.

Item	Description	Unit
8527929	Stainless steel	4 pcs



Dreve Otoplastik Cutter

Standard 5 1/2" stainless steel alligator forceps with serrated jaw. Not intended for invasive medical procedures.

Item	Description	Size	Unit
8530916	Pear-shaped	5.0 mm	1 pc
8530931	Conically-shaped	6.0 mm	1 pc



Dreve Otoplastik Ball Burr

Max. rotation speed (min-1) 30,000

Item	Size	Unit
8530914	Ø 2.1 mm	1 pc
8530915	Ø 3.5 mm	1 pc



Dreve Otoplastik Twist Drill

Max. rotation speed (min-1) 30,000.

Max. rotation speed (min-1) 30,001.

Max. rotation speed (min-1) 30,002.

Max. rotation speed (min-1) 30,003.

Item	Size	Unit
8530927	Ø 0.8 mm	1 pc
8530928	Ø 1.2 mm	1 pc
8530929	Ø 1.6 mm	1 pc
8530930	Ø 2.85 mm	1 pc



Westone Alligator Forceps

Standard 5 1/2" stainless steel alligator forceps with serrated jaw. Not intended for invasive medical procedures.

Item	Unit
8527901	1 pc

Westone Scissors, 3.5"

Sharp scissors suitable for general use in the clinic.

Part no.	Unit
8527905	1 pc



Westone Mirror Go Round

Mirror Go Round allows patients to easily view their ears to see how their hearing instruments should look when properly placed in the ear. Folds for easy transport and storage.

Item	Unit
8536591	1 pc



Dreve Otoplastik Demo Ear

Mockup made of flexible silicone. Ideal for demonstration purposes.

Item	Description	Unit
8536601	Left	1 pc
8536602	Right	1 pc



Dreve Otoplastik Silicone Glue SK 87

Item	Size	Unit
8532585	18 ml	1 pc



Dreve Otoplastik Superglue

Viscous

Item	Size	Unit
8530937	10 ml	1 pc



Premo Loctite 406 Instant Adhesive

An instant adhesive designed especially for the fast bonding of plastics and elastomeric materials where very fast fixturing is required. This product offers exceptionally low viscosity.

Item	Unit
8005067	1 pc

Hearing Aid Maintenance and Care

(HAS)



Westone Preformed Tubes

Item	Description	Unit
8527852	Standard #12	25 pcs
8527872	Medium #13	25 pcs
8527853	Thick #13	25 pcs
8527854	Super Thick #13	25 pcs

(HAS)



Westone Performed Disappear Tubing

Item	Description	Unit
8527855	Standard #12	10 pcs
8527856	Medium #13	25 pcs
8527857	Thick #13	25 pcs

(HAS)



Westone Tubing Retention System

Item	Description	Unit
8527927	Thick #13	25 pcs
8527928	Super Thick #13	25 pcs

(HAS)



Westone Tube-Lock With Tube

Item	Description	Unit
8527860	Thick #13	10 pcs
8527861	Super Thick #13	10 pcs

(HAS)



Westone Stay Dry-Tube

Item	Description	Unit
8527873	Medium #13	10 pcs
8527874	Thick #13	10 pcs

(HAS)



Westone Horn Tube

Item	Description	Size	Unit
8527934	Medium	3 mm	6 pcs
8527930	Thick	3 mm	6 pcs
8527931	Medium	4 mm	6 pcs

(HAS)



Westone Horn Tubes for TRS

Item	Size	Unit
8527932	3 mm	6 pcs
8527933	4 mm	6 pcs

(HAS)



Dreve Otoplastik Pre-bent Tubing

Item	Size	Unit
8530962	Ø2×3 mm.	10 pc

(HAS)



Westone Tube-lock Without Tube

Brass tube-lock for use with 13 thick and 13 super thick tubing.

Item	Description	Unit
8527858	Thick	1 pc
8527859	Super Thick	1 pc

(HAS)



Dreve Otoplastik Arc Tubing Connection

With retentions.

Item	Unit
8530938	10 pcs

(HAS)



Dreve Otoplastik Tubing

Item	Description	Unit
8530939	Ø2×3 mm	10 m

(HAS)



Dreve Otoplastik OtoVita Tubing Floss

For cleaning of micro tubing of open BTE hearing systems.

Item	Description	Unit
8530957	Micro Tubing	30 pcs

For cleaning of standard tubing (2×3 mm) of SE earmolds.

Item	Description	Unit
8530958	Tubing	30 pcs

Hearing Aid Maintenance and Care



Westone Excel Soft Grip Razor Knife

Great for trimming a variety of materials features a rear locking device and safety cap.

Item	Description	Unit
8536589	3/8" x 5"	1 pc



Westone Precision Screwdriver Set

Screwdriver set for making adjustments to hearing instruments or other precision work.

Item	Description	Unit
8527896	Four flat-tip screwdrivers and two Phillips tips.	6 pcs



Westone Precision Screwdriver Set

Screwdriver set for making adjustments to hearing instruments or other precision work.

Item	Description	Unit
8527897	Six flat-tip screwdrivers, three Phillips tips, one needle-point tip, and one magnet tip.	11 pcs



Surgical Blades and Handles

Suitable for cutting HA tubing.

Item	Description	Unit
8515427	Handle No. 3	1 pc
8515428	Blades No. 11	100 pcs



Westone Tubing Expander Plier

Stainless steel pliers with two-way lever spring allow you to quickly and easily stretch acoustic tubing before fitting it over earhooks, elbows, etc.

Item	Description	Unit
8527875	Curved	1 pc
8527876	Straight	1 pc



Westone Tweezers

For a variety of clinical uses. Good for working on hearing instruments and other small parts.


Item	Description	Unit
8527902	Curved Point Precision 6"	1 pc
8527903	Fine Tipped 4.75"	1 pc
8527904	Anti-Magnetic 4.25"	1 pc
8536603	Bayonet	1 pc



Westone End Nipper Pliers

A perfect tool for trimming tubing from the end of a canal after the earmold has been retubed.

Item	Unit
8527900	1 pc



Westone Needle Nose Pliers 5"

Needle nose pliers for a multitude of uses in the hearing healthcare clinic.

Item	Unit
8527898	1 pc



Westone Sprue Cutter, 4.5"

A unique design easily trims tubing flush to the earpiece.

Item	Unit
8527935	1 pc

Take care of your hearing aids

Take good care of your hearing aid. Your Hearing instrument is a high tech device needing regular care in order to provide you the best performance every day.



Audinell Cleansing Tablets

Effervescent tablets are specially designed to clean earmolds. The soaking solution contains a detergent and active oxygen to dissolve earwax while maintaining the effectiveness of the hearing aids.

Item	Unit
8520148	Pack/4 pcs



Audinell Dehydrating Capsules

Drying tablets are used for the dehumidification of hearing aids and earmolds, preventing any issues related to humidity and ensuring the long-term effectiveness of hearing aids. These translucent capsules are packets and contain highly absorbent orange silica gel.

Item	Unit
8520149	Pack/4 pcs



Audinell Drying Cup

Individual Audinell hearing aid drying cup only for use with drying tables.

Item	Unit
8520150	1 pc



Audinell Cleaning Cup

The cleaning cup is an essential accessory to be used with cleaning tablets. It provides an adequate environment to ensure efficient cleaning of the earmolds through the action of the cleaning tablets.

Item	Unit
8520151	1 pc



Audinell OtoFloss

Special Cleansing threads for sound tubes.

Item	Unit
8520152	100 pcs



Audinell PerfectDry Lux

The PerfectDry Lux is a Medical Device for the maintenance and storage of all kind of hearing aids. It's modern and sleek design combined with the best technology makes it exceptionally convenient to use. It's efficiency in dehumidifying in only 30 minutes has been proved by many independent tests. Height: 21 x Width: 70 x Depp: 51 mm.

Item	Unit
8520154	1 pc

Hearing Aid Maintenance and Care



Audinell Cleaning Spray

The cleaning spray dissolves earwax while maintaining the effectiveness of hearing aids. It comes in a plastic bottle with a valve, a nozzle and a brush. It has an anti-rotation spray button for increased comfort when in use.

Item	Unit
8520147	1 pc



Audinell Natural Lubricant

The natural lubricant is recommended for those with newly fitted hearing aids or for everyday comfort. Made using sweet almond oil, it makes it easier to insert the hearing aid in the ear canal. It also works to prevent irritation caused by the daily wearing of hearing aids.

Item	Unit
8520052	Bottle/1 pc



Audinell Drying Spray

The high-pressure spray is designed to ensure quick, effective drying and remove all dust from the hearing aid. The spray is dry and neutral and won't damage any components or acrylic resin.

Item	Unit
8520153	1 pc



Audinell Ear Spray

Audinell ear spray is a gentle isotonic washing solution containing surfactants. It doesn't contain propellant gas and won't cause any unpleasant cold sensations. The formula's chamomile extract soothes irritation caused by the daily wearing of hearing aids.

Item	Unit
8520063	Spray/1 pc



Audinell Dermoprotective Gel

The dermoprotective gel hydrates and lubricates the ear canal to make it easier to insert hearing aids, without damaging them, and improve tolerance. Made with a natural and moisturizing active ingredient known for its powerful anti-irritation properties, it soothes the ear canal without irritating the skin. Its skin compatibility has been tested by a clinical institute.

Item	Unit
8520053	Tube /1 pc



Audinell Earmold Air Blowers

Removes moisture and foreign matter from molds and tubing after cleaning to provide a clear opening from hearing aid to eardrum.


Item	Unit
8515562	1 pc



Audinell Cleaning/Disinfecting Wipes

Cleaning and disinfecting wipes, specially designed for surfaces, instruments, healthcare seats and medical devices. Microbiological properties: The solution's specific formulation has virucide actions (according to RKI, BGVV and DVV), bactericide actions (as per the EN 13727 standard) and fungicide actions (as per the EN 14562 standard) for a contact time of 1 minute.

Item	Description	Unit
8520146	Cleaning Wipes	30 pcs
8525161	Disinfecting Wipes	90 pcs



Impression Materials and Supplies

With impression materials, you can create customized molds of the ear canal, ensuring a comfortable and precise fit for the hearing aid. Our range of materials, such as silicone or foam, are tailored to ear anatomies and hearing aid specifications. Our supplies also include syringes, trays, and mixing tips. Together, these supplies play a crucial role in ensuring that your patient's hearing aid meets their expectations for comfort, performance, and satisfaction.

Impression Materials and Supplies



(HAS)

Dreve Otoplastik Mixing Tips For Double Cartridge.

Item	Description	Color	Unit
8522432	5,4 mm, 50ml	Pink	100 pcs
8522430	5,4 mm, 48 ml		100 pcs
8522431	6,3 mm, 48 ml		100 pcs
8522183	6,5 mm, 50 ml	Green	100 pcs



(HAS)

Dreve Otoplastik Injector Manual

Manual hand unit for dispensing all standard 48 ml and 50 ml double cartridges. The injector allows safe and clean application of impression materials with even pressure. Little force is required to inject the silicone.

Item	Description	Unit
8522434	For 50 ml	1 pc
8522108	For 48 ml	1 pc



(HAS)

Westone Foam PROS®

These revolutionary oto-dams are manufactured with a ventilation tube that relieves the air pressure that can build up in the ear canal while impression material is expressed into the ear, or when the cured impression is removed. OtoBlock - Slim PROS pressure relief Oto-dam system.

Item	Description	Unit
8536582	9 mm - Black	24 pcs
8536583	11 mm - Beige	24 pcs
8536584	13 mm - Blue	24 pcs



(HAS)

Westone Bite Blocks

Open jaw/bite block impressions for hearing protection devices.

Item	Color	Unit
8527850	White	25 pcs
8527849	White	100 pcs



(HAS)

Dreve Otoplastik Mixing Tray With dosing scale.

Item	Unit
8522105	1 pc



(HAS)

Dreve Otoplastik Mixing Pad

Item	Unit
8522106	100 sheets



(HAS)

Dreve Otoplastik Large Mixing Spatula Made of metal with plastic handle.

Item	Description	Unit
8522104	Metal	1 pc

Impression Materials and Supplies

HAS

Dreve Otoplastik OtoForm

This impression material features a pasty consistency while maintaining absolute thixotropic stability, ensuring it does not flow out of the ear. After curing, the impression can be easily removed, demonstrating high tensile strength. It is particularly suited for CIC impressions, very small auditory canals, and individuals sensitive to pain. The material is versatile, supporting both pressure-free and slightly pressure-forming impression techniques, making it ideal for fast, secure, and hygienic impression taking.

Item	Description	Size	Unit
8522433	A Flex	8x50 ml	1 pc
8536593	A Soft	8x50 ml	1 pc
8536594	A SoftX	8x50 ml	1 pc
8530889	Ak	8x50 ml	1 pc



HAS

Dreve Otoplastik FotoTec Repair

FotoTec® Repair is a light-curing one component resin in premium quality for the repair of ITE shells and BTE earmolds.

Item	Description	Unit
8536596	Blue transparent	1 pc
8536597	Transparent	1 pc
8536598	Red transparent	1 pc
8536599	Beige opaque 1	1 pc
8536595	Beige opaque 2	1 pc



HAS

Oticon A-ZOFT Impression Material

Premium quality and smooth initial consistency. Extremely soft, easy to remove and allows highest impression precision.

Item	Description
8515424	2x800 gram jars
8515426	40 foil bags (sets 2x8 g)



HAS

Dreve Otoplastik Otoform AK

Impression material, non-sticky and smooth with a slightly pressure-forming consistency. An even and controlled pressure is applied on the ear tissue. Otoform® Ak is advantageous for earmold fitting and sealing for less resilient ear tissue. Tubes (A+B), inc, dosing spoons.

Item	Size	Unit
8528835	2x544 ml	1 pc



HAS

Oticon Impression Syringe

Easy to fill, inject and clean. Provides a stable and easy injection.

Item	Size	Unit
8515445	4 mm	1 pc



HAS

Dreve Otoplastik Fotoplast Lack 3

For lacquering resp. "glazing" of all hard earmolds instead of mechanical polishing. Curing is done within 60 seconds (instant lacquer). Brilliant, tissue-friendly surface sealing. Biocompatible.

Item	Size	Unit
8530893	50 ml	1 pc



HAS

Dreve Otoplastik Fotoplast Gel

Light curing repair resin in thixotrope-stable consistency. Packed in tubes or dosing syringes with 3 long, flexible dosing tips. Working on ventings and repairing is very comfortable with the flexible dosing syringes.

Item	Description	Size	Unit
8530895	Clear transparent	20 ml	1 pc
8530896	Opaque H	20 ml	1 pc




Impression Materials and Supplies



(HAS)

Foam Impression Pads Slim

Item	Description	Color	Unit
8527869	Approx. 7 mm	White	50 pcs
8515429	Approx. 9.5 mm	Mint	50 pcs
8515430	Approx. 11 mm	White	50 pcs
8515423	Approx. 12.5 mm	Blue	50 pcs




(HAS)

Egger Impression Pads

Impression pads made of elastic foam with securing thread.

Item	Description	Color	Unit
8515560	12 mm	Pink	100 pcs




(HAS)

Dreve Otoplastik Impression Pads

Foam cone with safety thread to protect the ear drum during impression taking. The pads consist of flexible foamed material and seal the auditory canal by a safe fitting. Insertion with jaw tweezers and ear light. The impression pads CIC are suitable for children or narrow auditory canals.

Item	Description	Color	Unit
8522182	7 mm	Yellow	100 pcs
8515559	11 mm	Rose	100 pcs



(HAS)

Foam Impression Pads Slim

With white safety thread.

Item	Description	Color	Unit
8527870	Mini	White	50 pcs
8527871	Small	White	50 pcs
8515572	Medium	White	50 pcs
8515571	Small	Blue	50 pcs





Earplugs and Batteries

Send your patients home with ear plugs they can use for flying or swimming and make sure their hearing instrument is powered up with our selection of hearing instrument batteries.

Earplugs and batteries



Sanohra Earplugs

Sanohra Fly or Swim earplugs, suitable for children and adults.

Item	Description	Unit
8519351	Kids - Fly	1 set
8519352	Adult - Fly	1 set
8519353	Kids - Swim	1 set
8520143	Adult - Swim	1 set



(HAS)

Rayovac Hearing Instrument Batteries

Fits all types of Hearing Instruments.

Item	Description	Unit
8520134	Rayovac Extra Zap 10	6 pcs
8520135	Rayovac Extra Zap 13	6 pcs
8520136	Rayovac Extra Zap 312	6 pcs
8520137	Rayovac Extra Zap 675	6 pcs



(HAS)

Rayovac Battery Accessories

Item	Description	Unit
8515452	Battery tester	1 pc
8520138	Battery "carrier"	1 pc



Balance Supplies

Our range of balance supplies is designed to assist in conducting detailed analyses of balance function. Additionally, we provide a selection of consumables and spare parts to ensure that your equipment is always in top condition.

Balance Supplies



Difra Airstar Disposable specula

Item	Description	Color	Unit
8011457	2.5 mm spec.	Green	100 pcs
8011456	4 mm spec.	Blue	100 pcs




Sanibel Disposable VNG Goggle Foam Pads

Designed for:

Interacoustics VNG mask.

Item	Color	Unit
8011527	Black	24 pcs



Sanibel Disposable Ear Specula

Designed for:

Interacoustics AirFx™.

Item	Description	Color	Unit
8013013	2.75 mm	Black	34 pcs
8013014	4.25 mm	Black	34 pcs



Sanibel Water Irrigator Tips

Designed for:

Interacoustics Aqua Stim™.

Item	Color	Unit
8013018	Transparent	50 pcs



Infection Control

Preventing the spread of bacteria, dirt, and dust is essential for hygienic hearing testing.

Accessories such as gloves, shoe covers, and aprons play a critical role in maintaining high professional standards and creating a safe clinical environment during patient interactions. Explore more on the following pages to learn about our infection control offerings.

Infection Control



Nitrile Examination Gloves

To ensure total comfort, the 'thermo-elastic' material adapts to the hand's shape. Do not interact with petroleum oils, grease, solvents.

Item	Size	Color	Unit
8515551	Small	Light blue	100 pcs
8515552	Medium	Light blue	100 pcs
8515553	Large	Light blue	100 pcs
8515554	Extra large	Light blue	90 pcs



Disposable apron

Keep a professional appearance also after servicing a hearing instrument.

Item	Unit
8519348	100 pcs



Baaske Medical Keyboard Protective Foil

Cover for a standard PC keyboard, effectively protects against bacteria and germs. Service life: up to 3 months (depending on the application).

Item	Description	Unit
8500460	400 to 500 mm	1 pc



Dreve OtoVita Disinfection Spray

For disinfecting ear impressions and work surfaces to ensure a hygienic workplace.

Item	Description	Unit
8536605	750 ml	1 pc



Dreve Otovita Hand Gel

For hygienic hand washing according to DIN 1499 without water. The hand gel is characterized by a smooth consistency, short contact time and natural odor.


Item	Description	Unit
8536606	1l	1 pc



Dreve OtoVita Professional Disinfection Concentrate

For the disinfection of earmolds and instruments in the ultrasonic cleaning unit. We recommend to change the ultrasonic bath daily to assure an ideal disinfection.

Item	Description	Unit
8530888	1l	1 pc



Thermal Printer Paper and Labels

We offer a wide range of paper and labels for printers, specifically designed for use in hearing tests. Our products come in various sizes and formats to meet all your needs. Whether it's printing test results or organizing patient information, our high-quality products ensure clear prints and reliable performance.

Thermal Printer Paper and Labels



Thermal Paper
For **ERO-SCAN™/OtoRead™**
Width: 57 mm.
Inner diameter: 12 mm.
Outer diameter: 45 mm.

Item	Unit
8500200	1 roll/25 m



Thermal Paper LPQ80
For **GSI TympStar™ Pro.**
Width: 80 mm.
Inner diameter: 13 mm.
Outer diameter: 80 mm

Item	Unit
8509893	1 roll/80 m



API300 Thermal Paper
For **Titan.**
Width: 57 mm.
Inner diameter: 12 mm.
Outer diameter: 35 mm.

Item	Unit
8500230	1 roll/10 m




TPR10™ Thermal Paper
for **MTP10.**
Width: 79 mm.
Inner diameter: 12 mm.
Outer diameter: 45 mm

Item	Unit
8500220	1 roll/25 m




PRT MPT-II Printer
For **OtoRead.**
Width: 57 mm
Inner diameter: 8 mm
Outer diameter: 30 mm

Item	Unit
8013009	1 roll/10 m



TPR26™ Thermal Paper
For **Interacoustics.** Fits most popular audiometers and tympanometers.
Width: 112 mm.
Inner diameter: 12 mm.
Outer diameter: 45 mm.

Item	Unit
8500210	1 roll/25 m




Label Paper Roll HM-E200
Width: 56 mm,
Length: 60 mm.
R angle: 2.
Roller size:
Paper width: 57 mm.
Label spacing: 6 mm.
Core diameter: 13 mm.
Roller diameter: less or equal to 40 mm.

Item	Unit
8519836	1 roll/120 labels



Adhesive Thermal Paper LPQ80
For **GSI TympStar Pro.**
Width: 82 mm.
Inner diameter: 25 mm.
Outer diameter: 88 mm.

Item	Unit
8509894	1 roll/195 labels



MPT-III Thermal Paper
For **Sanibel MPT-III** printer which is supplied with **Interacoustics AT235.**
Width: 78 mm.
Inner diameter: 12 mm.
Outer diameter: 40 mm.








Item	Unit
8500610	1 roll/20 m



Thermal Paper
For the **Sanibel MPT-II** printer which is supplied with **MAICO ERO-SCAN™, GSI Corti™ and Interacoustics OtoRead™.**
Width: 57 mm.
Inner diameter: 9 mm.
Outer diameter: 40 mm.

Item	Unit
8029305	1 roll/20 m

List of Abbreviations

Symbol/ marking	Reference	Title	Description
	Not a standardized symbol	Auditory Brainstem Response	Auditory Brainstem Response (ABR) is a test that measures the electrical activity in the brainstem in response to auditory stimuli, commonly used to assess hearing, especially in infants or non-cooperative individuals.
	Not a standardized symbol	Audiometry	Audiometry is a method of assessing a person's hearing by measuring their ability to perceive sound, typically involving presenting various sounds at varying intensities and frequencies to determine a person's hearing threshold, or the lowest volume they can hear at different frequencies.
	Not a standardized symbol	Balance	Balance refers to the ability to maintain bodily stability by integrating sensory information from the inner ear with other senses.
	Not a standardized symbol	Hearing Aid Supplies	Hearing Aid Supplies refers to various materials, equipment, or accessories used in conjunction with hearing aids.
	Not a standardized symbol	Hearing Aid Fittings	Hearing Aid Fittings refers to the process of customizing a hearing aid to fit an individual user's ear. This often involves taking measurements of the ear and its canals to ensure a comfortable and effective fit.
	Not a standardized symbol	Impedance	Impedance measurements assess the condition of the middle ear.
	Not a standardized symbol	Otoacoustics Emissions	The OAE test assesses the integrity and function of the outer hair cells within the cochlea, and can be used for differential diagnosis of hearing conditions.



Sanibel Supply A/S • info@sanibelsupply.com • www.sanibelsupply.com



Disclaimer

This catalogue has been published for worldwide distribution. Some parts described may NOT BE AVAILABLE in your area. Your local distributor will be pleased to give you details. For advice on compatibility and instructions of use, please consult your operation manual or local distributor. Sanibel Supply A/S assumes no legal liability or responsibility for the accuracy, typing errors or mistakes, correctness or completeness of any information in this flyer. The information in this flyer was correct to the best of our knowledge at the time of printing. Sanibel Supply A/S reserves the right to make changes at any time in specification of products described in this flyer without notice or obligation. MADSEN AccuScreen® and MADSEN OTOflex® are registered trademarks of Otometrics A/S. Echo Screen®, AuDx®, Navigator® Pro, and ABaer® are registered trademarks of Natus Medical Incorporated. Sanibel® is a registered trademark of Sanibel Supply®. All other instrument brands or model names are copyrighted or trademarked by their respective owners and are used only for reference.